

DISCRETE PYROS

MONITORS

ENERGY DETECTORS

POWER DETECTORS

HIGH POWER SOLUTIONS


PHOTO DETECTORS

THZ DETECTORS

OEM DETECTORS

SPECIAL PRODUCTS

BEAM DIAGNOSTICS



KEY FEATURES

- BROAD SPECTRAL RESPONSE**
From 0.1 to 1000 μm
- EASY TO INTEGRATE FORMAT**
T05 and T08 packages make the QS detectors small and easy to integrate in an existing system
- LARGE AREA SENSORS**
5 mm \varnothing and 9 mm \varnothing diameter pyroelectric sensors make optical alignment easier
- SEVERAL IR WINDOWS IN OPTION**
 - Quartz: 0.2 – 3.0 μm
 - Barium Fluoride: 0.2 – 17.5 μm
 - Sapphire: 0.1 – 7.0 μm
 - Silicon: 1.1 – 9.0 μm and 50 – 1000 μm
 - AR Germanium: 8 – 14 μm

AVAILABLE MODELS

5 families of products to choose from:

- QS-L Discrete Pyro Detectors, Low Noise Level
- QS-H Discrete Pyro Detectors, High Average Power
- QS-VL Hybrid Pyro Detectors, Voltage Mode, Low Noise Level
- QS-IF Hybrid Pyro Detectors, Current Mode, Fast Response
- QS-IL Hybrid Pyro Detectors, Current Mode, Low Noise Level

ACCESSORIES



QS-I-TEST

Evaluation Test Box (current)



QS-V-TEST

Evaluation Test Box (voltage)



Permanent IR Windows
(Various types available)



Pelican Carrying Case

SEE ALSO

TECHNICAL DRAWINGS	134
ABSORPTION CURVES	136
LIST OF ALL ACCESSORIES	188
APPLICATION NOTES	
COMPENSATING CURRENT MODE AMPLIFICATION USING QS-I-TEST	201925
HOW TO HANDLE SENSITIVE PYROELECTRIC DETECTORS	202181
THERMAL SATURATION IN HYBRID PYROELECTRIC DETECTORS	201926
HOW THEY WORK: QS-I-TEST & QS-V-TEST	201927
HIGHEST PERFORMANCE WITH QS DETECTORS	201928
QS DETECTORS PIN-OUTS & DESCRIPTIONS	201931
QS-I-TEST SPECIFICATIONS	202187

QS-L



Approved or in the process of being approved*

Discrete Pyro Detectors, Low Noise Level

SPECIFICATIONS

	QS1-L	QS2-L	QS3-L	QS5-L	QS9-L
CURRENT RESPONSIVITY	1 $\mu\text{A/W}$	0.5 $\mu\text{A/W}$	0.5 $\mu\text{A/W}$	0.25 $\mu\text{A/W}$	0.25 $\mu\text{A/W}$
EFFECTIVE APERTURE	1 mm \emptyset	2 mm \emptyset	3 mm \emptyset	5 mm \emptyset	9 mm \emptyset
PACKAGE	T05	T05	T05	T05	T08
MEASUREMENT CAPABILITY					
Spectral Range	0.1 - 1000 μm	0.1 - 1000 μm	0.1 - 1000 μm	0.1 - 1000 μm	0.1 - 1000 μm
Max Average Power	50 mW	50 mW	50 mW	50 mW	50 mW
Capacitance (at 1000 Hz)	15 pF	22 pF	60 pF	90 pF	250 pF
Current Responsivity (at 630 nm)	1 $\mu\text{A/W}$	0.5 $\mu\text{A/W}$	0.5 $\mu\text{A/W}$	0.25 $\mu\text{A/W}$	0.25 $\mu\text{A/W}$
Thermal Frequency (3 dB)	3.5 Hz	1.6 Hz	0.8 Hz	0.5 Hz	0.25 Hz
Temperature Coefficient	0.2%/°C	0.2%/°C	0.2%/°C	0.2%/°C	0.2%/°C
PHYSICAL CHARACTERISTICS					
Effective Aperture	1 mm \emptyset	2 mm \emptyset	3 mm \emptyset	5 mm \emptyset	9 mm \emptyset
Package	T05	T05	T05	T05	T08
Sensor	Pyroelectric	Pyroelectric	Pyroelectric	Pyroelectric	Pyroelectric
Absorber	MT	MT	MT	MT	MT
Dimensions	8.3 \emptyset x 6.4D mm	8.3 \emptyset x 6.4D mm	8.3 \emptyset x 6.4D mm	8.3 \emptyset x 6.4D mm	13.6 \emptyset x 6.4D mm
Weight	1.0 g	1.0 g	1.0 g	1.0 g	1.5 g
ORDERING INFORMATION					
Product Name	QS1-L	QS2-L	QS3-L	QS5-L	QS9-L
Product Number	201657	201659	201662	201664	201666

Specifications are subject to change without notice

QS-H



Approved or in the process of being approved*

Discrete Pyro Detectors, High Average Power

SPECIFICATIONS

	QS1-H	QS2-H	QS3-H	QS5-H	QS9-H
CURRENT RESPONSIVITY	500 mW	500 mW	500 mW	500 mW	500 mW
EFFECTIVE APERTURE	1 mm Ø	2 mm Ø	3 mm Ø	5 mm Ø	9 mm Ø
PACKAGE	T05	T05	T05	T05	T08
MEASUREMENT CAPABILITY					
Spectral Range	0.1 - 1000 µm	0.1 - 1000 µm	0.1 - 1000 µm	0.1 - 1000 µm	0.1 - 1000 µm
Max Average Power	500 mW	500 mW	500 mW	500 mW	500 mW
Capacitance (at 1000 Hz)	3 pF	12 pF	30 pF	90 pF	250 pF
Current Responsivity (at 630 nm)	0.25 µA/W	0.25 µA/W	0.25 µA/W	0.25 µA/W	0.25 µA/W
Thermal Frequency (3 dB)	5 Hz	5 Hz	5 Hz	5 Hz	5 Hz
Temperature Coefficient	0.2%/°C	0.2%/°C	0.2%/°C	0.2%/°C	0.2%/°C
PHYSICAL CHARACTERISTICS					
Effective Aperture	1 mm Ø	2 mm Ø	3 mm Ø	5 mm Ø	9 mm Ø
Package	T05	T05	T05	T05	T08
Sensor	Pyroelectric	Pyroelectric	Pyroelectric	Pyroelectric	Pyroelectric
Absorber	MT	MT	MT	MT	MT
Dimensions	8.3Ø x 6.4D mm	8.3Ø x 6.4D mm	8.3Ø x 6.4D mm	8.3Ø x 6.4D mm	13.6Ø x 6.4D mm
Weight	1.0 g	1.0 g	1.0 g	1.0 g	1.5 g
ORDERING INFORMATION					
Product Name	QS1-H	QS2-H	QS3-H	QS5-H	QS9-H
Product Number	201658	201661	201663	201665	201667

Specifications are subject to change without notice

QS-VL

Hybrid Pyro Detectors, Voltage Mode, Low Noise Level



Approved or in the process of being approved*

SPECIFICATIONS

	QS1-VL	QS2-VL	QS3-VL	QS5-VL	QS9-VL
VOLTAGE RESPONSIVITY	900 V/W	200 V/W	90 V/W	25 V/W	15 V/W
CURRENT RESPONSIVITY	1 μ A/W	0.5 μ A/W	0.5 μ A/W	0.25 μ A/W	0.25 μ A/W
EFFECTIVE APERTURE	1 mm \emptyset	2 mm \emptyset	3 mm \emptyset	5 mm \emptyset	9 mm \emptyset
PACKAGE	T05	T05	T05	T05	T08

MEASUREMENT CAPABILITY

Spectral Range	0.1 - 1000 μ m	0.1 - 1000 μ m	0.1 - 1000 μ m	0.1 - 1000 μ m	0.1 - 1000 μ m
Max Average Power	50 mW	50 mW	50 mW	50 mW	50 mW
Noise Equivalent Power ^a	3×10^{-10} W/(Hz) ^{1/2}	6×10^{-10} W/(Hz) ^{1/2}	1×10^{-9} W/(Hz) ^{1/2}	2×10^{-9} W/(Hz) ^{1/2}	5×10^{-9} W/(Hz) ^{1/2}
Detectivity ^a	2.9×10^8 cm(Hz) ^{1/2} /W	3.0×10^8 cm(Hz) ^{1/2} /W	2.7×10^8 cm(Hz) ^{1/2} /W	2.2×10^8 cm(Hz) ^{1/2} /W	1.6×10^8 cm(Hz) ^{1/2} /W
Capacitance (at 1000 Hz)	15 pF	22 pF	60 pF	90 pF	250 pF
Current Responsivity (at 630 nm)	1 μ A/W	0.5 μ A/W	0.5 μ A/W	0.25 μ A/W	0.25 μ A/W
Voltage Responsivity ^b	900 V/W	200 V/W	90 V/W	25 V/W	15 V/W
Thermal Frequency (3 dB)	3.5 Hz	1.6 Hz	0.8 Hz	0.5 Hz	2.5 Hz
Load Resistor	300 G Ω	300 G Ω	100 G Ω	100 G Ω	100 G Ω
Supply Voltage	+9 to +15 V	+9 to +15 V	+9 to +15 V	+9 to +15 V	+9 to +15 V

PHYSICAL CHARACTERISTICS

Effective Aperture	1 mm \emptyset	2 mm \emptyset	3 mm \emptyset	5 mm \emptyset	9 mm \emptyset
Package	T05	T05	T05	T05	T08
Sensor	Pyroelectric	Pyroelectric	Pyroelectric	Pyroelectric	Pyroelectric
Absorber	MT	MT	MT	MT	MT
Dimensions	8.3 \emptyset x 6.4D mm	8.3 \emptyset x 6.4D mm	8.3 \emptyset x 6.4D mm	8.3 \emptyset x 6.4D mm	13.6 \emptyset x 6.4D mm
Weight	1.0 g	1.0 g	1.0 g	1.0 g	1.5 g

ORDERING INFORMATION

Product Name	QS1-VL	QS2-VL	QS3-VL	QS5-VL	QS9-VL
Product Number	201673	201674	201676	201677	201678

Specifications are subject to change without notice

a. 630 nm, 5 Hz, 1 Hz Bandwidth

b. 630 nm, 15 Hz



QS-V-TEST EVALUATION TEST BOX

QS-V-TEST	
Batteries	+9V
R _i Resistors	10 ⁵ - 10 ¹⁰ Ω
C _i Compensating	NO
Package	101.6H x 127W x 58.4D
Optical Mount	1/4-20 Threaded
Front Bezel	SM1 (1.035-40)
Product Number	201694

* For details, contact your Gentec-EO representative

QS-IF



Approved or in the process of being approved*

Hybrid Pyro Detectors, Current Mode, Fast Response

SPECIFICATIONS

	QS1-IF	QS2-IF	QS3-IF	QS5-IF	QS9-IF
VOLTAGE RESPONSIVITY	100 V/W	50 V/W	50 V/W	25 V/W	25 V/W
CURRENT RESPONSIVITY	1 μ A/W	0.5 μ A/W	0.5 μ A/W	0.25 μ A/W	0.25 μ A/W
EFFECTIVE APERTURE	1 mm \emptyset	2 mm \emptyset	3 mm \emptyset	5 mm \emptyset	9 mm \emptyset
PACKAGE	T05	T05	T05	T05	T08

MEASUREMENT CAPABILITY

Spectral Range	0.1 - 1000 μ m	0.1 - 1000 μ m	0.1 - 1000 μ m	0.1 - 1000 μ m	0.1 - 1000 μ m
Max Average Power	50 mW	50 mW	50 mW	50 mW	50 mW
Noise Equivalent Power ^a	5×10^{-8} W/(Hz) ^{1/2}	8×10^{-8} W/(Hz) ^{1/2}	8×10^{-8} W/(Hz) ^{1/2}	1.6×10^{-7} W/(Hz) ^{1/2}	1.6×10^{-7} W/(Hz) ^{1/2}
Detectivity ^a	1.8×10^6 cm(Hz) ^{1/2} /W	2.2×10^6 cm(Hz) ^{1/2} /W	3.3×10^6 cm(Hz) ^{1/2} /W	2.8×10^6 cm(Hz) ^{1/2} /W	5.0×10^6 cm(Hz) ^{1/2} /W
Capacitance (at 1000 Hz)	15 pF	22 pF	60 pF	90 pF	250 pF
Current Responsivity (at 630 nm)	1 μ A/W	0.5 μ A/W	0.5 μ A/W	0.25 μ A/W	0.25 μ A/W
Voltage Responsivity ^b	100 V/W	50 V/W	50 V/W	25 V/W	25 V/W
Thermal Frequency (3 dB)	3.5 Hz	1.6 Hz	0.8 Hz	0.5 Hz	0.25 Hz
Feedback Resistor	100 M Ω	100 M Ω	100 M Ω	100 M Ω	100 M Ω
Supply Voltage	\pm 12 V	\pm 12 V	\pm 12 V	\pm 12 V	\pm 12 V

PHYSICAL CHARACTERISTICS

Effective Aperture	1 mm \emptyset	2 mm \emptyset	3 mm \emptyset	5 mm \emptyset	9 mm \emptyset
Package	T05	T05	T05	T05	T08
Sensor	Pyroelectric	Pyroelectric	Pyroelectric	Pyroelectric	Pyroelectric
Absorber	MT	MT	MT	MT	MT
Dimensions	8.3 \emptyset x 6.4D mm	8.3 \emptyset x 6.4D mm	8.3 \emptyset x 6.4D mm	8.3 \emptyset x 6.4D mm	13.6 \emptyset x 6.4D mm
Weight	1.0 g	1.0 g	1.0 g	1.0 g	1.5 g

ORDERING INFORMATION

Product Name	QS1-IF	QS2-IF	QS3-IF	QS5-IF	QS9-IF
Product Number	201679	201680	201681	201682	201683

Specifications are subject to change without notice

a. 630 nm, 15 Hz, largeur de bande de 1 Hz

b. 630 nm, 15 Hz



QS-I-TEST EVALUATION TEST BOX

QS-I-TEST	
Batteries	+9V/-9V
R _i Resistors	10 ⁵ - 10 ¹⁰ Ω
C _i Compensating	Yes
Package	101.6H x 127L x 58.4P
Optical Mount	1/4-20 Threaded
Front Bezel	SM1 (1.035-40)
Product Number	201693

* For details, contact your Gentec-EO representative

QS-IL



Approved or in the process of being approved*

Hybrid Pyro Detectors, Current Mode, Low Noise Level

SPECIFICATIONS

	QS1-IL	QS2-IL	QS3-IL	QS5-IL	QS9-IL
VOLTAGE RESPONSIVITY	50 kV/W	25 kV/W	25 kV/W	13 kV/W	13 kV/W
CURRENT RESPONSIVITY	1 μ A/W	0.5 μ A/W	0.5 μ A/W	0.25 μ A/W	0.25 μ A/W
EFFECTIVE APERTURE	1 mm \emptyset	2 mm \emptyset	3 mm \emptyset	5 mm \emptyset	9 mm \emptyset
PACKAGE	T05	T05	T05	T05	T08

MEASUREMENT CAPABILITY					
Spectral Range	0.1 - 1000 μ m	0.1 - 1000 μ m	0.1 - 1000 μ m	0.1 - 1000 μ m	0.1 - 1000 μ m
Max Average Power	50 mW	50 mW	50 mW	50 mW	50 mW
Noise Equivalent Power ^a	8×10^{-10} W/(Hz) ^{1/2}	2×10^{-9} W/(Hz) ^{1/2}	2×10^{-9} W/(Hz) ^{1/2}	6×10^{-9} W/(Hz) ^{1/2}	6×10^{-9} W/(Hz) ^{1/2}
Detectivity ^a	1.1×10^8 cm(Hz) ^{1/2} /W	9.0×10^7 cm(Hz) ^{1/2} /W	1.3×10^8 cm(Hz) ^{1/2} /W	7.0×10^7 cm(Hz) ^{1/2} /W	1.3×10^8 cm(Hz) ^{1/2} /W
Capacitance (at 1000 Hz)	15 pF	22 pF	60 pF	90 pF	250 pF
Current Responsivity (at 630 nm)	1 μ A/W	0.5 μ A/W	0.5 μ A/W	0.25 μ A/W	0.25 μ A/W
Voltage Responsivity ^b	50 kV/W	25 kV/W	25 kV/W	13 kV/W	13 kV/W
Thermal Frequency (3 dB)	3.5 Hz	1.6 Hz	0.8 Hz	0.5 Hz	0.25 Hz
Feedback Resistor	100 G Ω	100 G Ω	100 G Ω	100 G Ω	100 G Ω
Supply Voltage	± 5 to ± 12 V	± 5 to ± 12 V	± 5 to ± 12 V	± 5 to ± 12 V	± 5 to ± 12 V

PHYSICAL CHARACTERISTICS					
Effective Aperture	1 mm \emptyset	2 mm \emptyset	3 mm \emptyset	5 mm \emptyset	9 mm \emptyset
Package	T05	T05	T05	T05	T08
Sensor	Pyroelectric	Pyroelectric	Pyroelectric	Pyroelectric	Pyroelectric
Absorber	MT	MT	MT	MT	MT
Dimensions	8.3 \emptyset x 6.4D mm	8.3 \emptyset x 6.4D mm	8.3 \emptyset x 6.4D mm	8.3 \emptyset x 6.4D mm	13.6 \emptyset x 6.4D mm
Weight	1.0 g	1.0 g	1.0 g	1.0 g	1.5 g

ORDERING INFORMATION					
Product Name	QS1-IL	QS2-IL	QS3-IL	QS5-IL	QS9-IL
Product Number	201684	201685	201686	201687	201688

Specifications are subject to change without notice

- a. 630 nm, 5 Hz, 1 Hz Bandwidth
b. 630 nm, 15 Hz




QS-I-TEST EVALUATION TEST BOX

QS-I-TEST	
Batteries	+9V/-9V
R, Resistors	10^5 - 10^{10} Ω
C, Compensating	Yes
Package	101.6H x 127L x 58.4P
Optical Mount	1/4-20 Threaded
Front Bezel	SM1 (1.035-40)
Product Number	201693


* For details, contact your Gentec-EO representative

DISCRETE PYROS



PYROELECTRIC THERMAL DETECTORS

Our pyroelectric detectors are a class of room temperature thermal detectors that produce a current output that is directly proportional to the rate of change of temperature when exposed to a source of radiation. They are best described by an AC current source, capacitor and resistor. Their current output is governed by the equation $I = p(T) \cdot A \cdot dT/dt$, where I is current, $p(T)$ is the Pyro Coefficient, A is the area as defined by the front electrode, and dT/dt is the rate of temperature change of the pyro crystal. The advantages of a pyroelectric detector over other IR detectors are: room temperature operation, broad spectral response, high sensitivity (D^*) and fast response (sub-nsec into 50Ω).




QS-L (left) and QS-H (right) Pin-Outs

QS-L AND QS-H DISCRETE PYROS

Our passive Discrete Pyroelectric Detectors range from 1 to 9 mm in diameter and are provided in two configurations: high sensitivity or high average power. They present a pyroelectric detector element covered with our metallic coating (MT) and are packaged in a miniature TO-5 or TO-8 can. The diagram shown left identifies the Pin-out for both types of detectors. Our organic black coating (BL), increases the optical absorption and helps flatten the spectral response. We also offer a number of permanent IR Windows that can be added to the TO can. These discrete pyro detectors are ideal for pulsed laser applications.

QS-VL VOLTAGE MODE HYBRID PYROS


Pyroelectric Detectors are high impedance ($>10^{13} \Omega$) devices that require use in an impedance converting circuit when trying to achieve the highest sensitivity (high D^*). Our QS-VL Series detectors include our pyroelectric element mated to an ultra-low noise FET in a source follower circuit. These are in turn packaged in a miniature TO-5 or TO-8 can. The equivalent circuit and Pin-out for this series are shown at the left. They are also available in sizes ranging from 1 to 9 mm diameter. These models are ideal for analytical instrumentation applications like Broadband IR Radiometers, Optical Pyrometer, and/or FTIR Spectrometers.



QS-VL Pin-Out


QS-IF AND QS-IL CURRENT MODE HYBRID PYROS

These detectors offer high gain ($>10^5$ V/W) and/or high bandwidth (>10 MHz). In this configuration, the pyroelectric detector element is combined to a low noise operational amplifier. The QS-IL models are designed for high performance at low to medium frequencies, while the QS-IF models offer good performance at medium to high frequencies. These detectors are very easy to use. Simply supply the ± 10 to 15 V to power the operational amplifier and add an external resistor, if required, to adjust the bandwidth and you are ready to measure pulsed, modulated or chopped sources, from nJ to mJ and nW to W. These detectors also make great candidates for any variety of broadband analytical instruments or laser measurement products.



QS-IF and QS-IL Pin-Out

DISCRETE PYROS




QS-VL and QS-IL Circuitry

VOLTAGE OUTPUT VS. FREQUENCY

Our QS-VL and QS-IL Hybrid Detectors are designed to maximize voltage output at low frequencies and therefore include load and feedback resistors in the 100 GΩ to 300 GΩ range. They are also designed into 8-pin TO packages that allow the addition of an “external resistor” to lower the output and increase the bandwidth. The circuit diagram at the left shows a typical hook up for our QS5-IL detector (with our MT coating), using external resistors and capacitors. Our QS-IF series, on the other hand, are designed for high bandwidth applications and therefore include a smaller feedback resistor of 100 MΩ. For expert help on designing a detector circuit please contact us info@gentec-eo.com.

OPERATION IN POWER MEASUREMENT MODE


When using our QS-IL Hybrid Detector to measure the Power (in Watts) of your CW or High Rep Rate source (Quasi-CW), you will need to employ an optical chopper. The diagram at the left shows the typical voltage output of a QS5-IL when used with our QS-I-TEST evaluation test box. Note that the voltage output is an approximate “square wave” whose rise and fall times are governed by the RC time constant of the circuit. The optical power is directly proportional to the peak voltage minus the baseline voltage. We calibrate these devices when operating in this mode.



Typical QS-IL Voltage Output in Power Measurement Mode

OPERATION IN ENERGY MEASUREMENT MODE


Our Pyroelectric Detectors are an ideal choice when measuring the performance of your pulsed laser in the range of nJ to mJ, across the full spectrum! The scope trace at the left represents the typical output from a QS9-IL, when used with our QS-I-TEST set up as an integrating Joulemeter. Note the fast rise to a peak and then slower decay governed by the RC time constant selected for the integrating circuit. In this configuration you can measure absolute pulse energy, rep rate, and pulse-to-pulse stability. The maximum pulse width of your source is determined by the RC time constant you select and there is no limit as to how short the pulse can be!



Typical QS-IL Voltage Output in Energy Measurement Mode

BROAD SPECTRAL RESPONSE

Unlike photoconductive and photovoltaic detectors, our Pyroelectric Thermal Detectors are not limited to a small part of the electromagnetic spectrum. They are truly broad spectrum detectors, sensitive from 0.1 μm to 3000 μm (EUV, FAR IR, and THz). Any and all radiation absorbed by our coatings or pyro crystal will result in a measurable signal. The two plots at the left show the relative spectral response of detectors with MT and BL coatings. Note that the well documented, NIST traceable calibrated portion of these curves runs from 0.25 μm to 15 μm. There are currently no traceable optical standards for measurements > 15 μm.



Absorption Curves of QS Pyroelectric Detectors